



Antibody Specification Sheet

Anti-human CD80 - Purified



Catalog reference: AS31-P
Size: 0.1 mg

Antibody Information:

Antigen: recombinant expressed extra-cellular domain of human CD80.

Ig Class: mIgG1.

Specificity: recognizes the extra-cellular domain of human CD80 (also called B7-1). Does not bind the extracellular domain of CD86.

Antibody Source: monoclonal antibody from BALB/c-derived hybridoma AS31

Production: *in vitro* cell culture.

Purification: Protein A agarose chromatography.

Purity: ≥95%.

Formulation: provided as a 0.2 µm sterile-filtered solution in phosphate buffered saline with 0.05% sodium azide.

Recommendations: Non-blocking antibody.

Storage conditions: store undiluted at 4°C..

Applications:

Neutralization: AS31 does not to block the interaction of soluble hCD80 with human CD152 (CTLA4).

Flow Cytometry:
(Please note - this application currently in development)

General Information:

Cell Distribution: CD80 is rapidly induced on the surface of activated B cells, T cells or monocytes and is expressed at high levels on dendritic cells.

Function: CD80 is one of two ligands (the other CD86) for CTLA4 and CD28. CD80 acts as co-stimulatory molecule in eliciting T-cell help during antigen presentation. Antigen presentation in the absence of sufficient costimulation involving CD86/CD80 can induce tolerance. CD80 appears to play a role distinct from CD86 in T helper cell differentiation.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Antibody Solutions

1206 Charleston Road, Mountain View, CA 94043, USA

Tel: 650-938-4300, FAX: 650-938-4390

E-mail: solns@antibody.com