



Antibody Specification Sheet

Anti-human CD152 - Purified



Catalog reference: AS32-P
Size: 0.1 mg

Antibody Information:

Antigen: recombinant expressed extra-cellular domain of human CD152.

Ig Class: mIgG1.

Specificity: recognizes the extra-cellular domain of human CD152 (also called CTLA-4).

Antibody Source: monoclonal antibody from BALB/c-derived hybridoma AS32.

Production: *in vitro* cell culture.

Purification: Protein A agarose chromatography.

Purity: ≥95%.

Formulation: provided as a 0.2 µm sterile-filtered solution in phosphate buffered saline

Recommendations: Neutralization.

Storage conditions: store undiluted at 4°C.

Applications:

Neutralization: AS32 has been shown to block the interaction of soluble hCD152 with hCD80.

General Information:

Cell Distribution: CD152 is expressed on most T lymphocytes. The level of expression is activation-dependent.

Function: CD152 is one of two ligands (the other CD28) for CD80 and CD86. CD152 acts as co-stimulatory molecule in eliciting T-cell help during antigen presentation.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Antibody Solutions

1206 Charleston Road, Mountain View, CA 94043, USA

Tel: 650-938-4300, FAX: 650-938-4390

E-mail: solns@antibody.com