



Antibody Specification Sheet

Anti-Human Ciliary NeuroTrophic Factor (hCNTF) - Serum

Catalog reference: PS2-S

Size: 0.1 ml



Antibody Information:

Antigen: human recombinant ciliary neurotrophic factor (hCNTF).

Ig Class: rb IgG.

Form: polyclonal serum.

Specificity: recognizes human CNTF. Does not

recognize a panel of unrelated human cytokines.

Antibody Source: from immunized rabbits.

Formulation with 0.05% sodium azide as a preservative.

Recommendations:

Storage conditions:

store undiluted at 4 °C or in aliquots at ≤ -20 °C.

Applications:

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Antibody Solutions

1130 Mountain View Alviso Road, Sunnyvale, CA 94089, USA

Tel: +1-650-938-4300, FAX: +1-408-734-0400

E-mail: solutions@antibody.com